

# Reducing Patient Identification Errors

while enhancing lean workflow

## Thermo Scientific PrintMate AS Cassette Printer

Eliminate waste and streamline your labeling process by printing cassettes on-demand with the Thermo Scientific™ PrintMate™ AS Cassette Printer at the grossing station.

- Prints each cassette in 10 seconds or less
- Small footprint allows the printer to fit easily next to the grossing station
- Print only the cassettes you need, right at the point of use
- Print barcodes on cassettes to automate downstream processes
- Eliminate manual handwriting processes



**New Product**



### Are specimen identification errors a concern in your laboratory?

The small footprint of the Thermo Scientific™ SlideMate™ AS Slide Printer brings the convenience of on-demand slide printing directly to the microtome. The SlideMate AS Slide Printer drastically reduces the chance for patient identification errors by allowing the laboratory to print accurate information at the moment the slide is required.

**Thermo**  
SCIENTIFIC



### Order Information

Description	Order Number
PrintMate AS 150 Cassette Printer	A83800001
PrintMate AS 450 Cassette Printer	A83800002
PrintMate AS 900 Cassette Printer	A83800003
Includes PrintMate AS, USB cable, power cable, 2, 6 or 12 cassette hoppers, one roll of tape (installed), LabWriter software CD, operator guide	
PrintMate AS 150 Cassette Printer	A83800004
PrintMate AS 450 Cassette Printer	A83800005
PrintMate AS 900 Cassette Printer	A83800006
Includes PrintMate AS, USB cable, power cable, 2, 6 or 12 cassette hoppers, one roll of tape (installed), Thermo Scientific™ LabWriter™ software CD, operator guide, all-in-one PC, digital imager	

### PrintMate AS Accessories

Description	Order Number
Collection system (includes cable, 4 trays for cassettes, output guides)	A79510014
Collection trays (box of 4)	A79510022
Hot foil tape, box of 5 rolls	A79510013
All-in-one touch screen PC	B81300033
Digital imager (linear/2D barcode reader)	B35310050
Cassette hopper (ea)	A79510024

### Specifications

#### Dimensions: (D x W x H):

PrintMate AS 900 w/o cassette tubes	15 x 17 x 11.4 in. (38 x 43.5 x 29 cm)
PrintMate AS 900 w/ cassette tubes	15 x 17 x 32 in. (38 x 43.5 x 81 cm)

#### Weight:

PrintMate AS 900 w/o cassette tubes	31 lb (14 kg)
PrintMate AS 900 w/ cassette tubes	37.5 lb (17 kg)

<b>Power Requirements:</b>	100–240 VAC, 50/60 Hz, 150 VA (512 BTU/h)
----------------------------	--

Print media type:	Thermal transfer ribbon
Print media capacity:	Approximately 20,000 cassettes per roll
Print Speed:	10 seconds per cassette
Connectivity:	USB
Print media color:	Black

#### Cassette Capacity:

PrintMate AS 150	150 cassettes
PrintMate AS 450	450 cassettes
PrintMate AS 900	900 cassettes

To learn more about how on-demand cassette printing can impact in your laboratory, call: (800) 522-7270

[www.thermoscientific.com/pathology](http://www.thermoscientific.com/pathology)

© 2014 Thermo Fisher Scientific Inc. All rights reserved. All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries.

#### Anatomical Pathology

4481 Campus Drive  
Kalamazoo, MI 49008  
UNITED STATES

+1 (800) 522-7270

Tudor Road, Manor Park  
Runcorn, WA7 1TA  
UNITED KINGDOM

+44 (0) 800 018 9396  
+44 (0) 1928 534 000